

<u>Equipment Item</u>	<u>Treatment</u>	<u>Comments</u>
Ozone Generator	All	Converts medical oxygen (O2) to Ozone (O3)
Oxygen Tank *	All	870 Medical grade with toggle valve (Size D recommended)
Oxygen Regulator	All	870 Medical pediatric/low-flow to fit an 870 toggle valve tank
Silicone Tubing	All	Required to connect the equipment together
Connectors (Luer/Push)	All	Required to connect the tubing to equipment
Oil Bubbler	Respiratory	Makes ozone safe for breathing
Syringe & Catheter	Rectal	Used to deliver rectal ozone to your pet
Insufflation Bag w/catheter	Rectal	Used to deliver rectal ozone to yourself
Limb Bag	Limb	Bags a limb for wound or skin treatment
Water Bubbler	Water	Creates ozonated water for drinking or disinfecting
3-way Stop Cock	Rectal	Easily switches ozone to either syringe/bag or destruct
Ozone Destruct	Rectal/Limb/Water	Prevents breathing of pure ozone during treatment
Check Valve/Air Trap	Respiratory/Water	Prevents oil/water from backing up into ozone generator